

TOURER (ST4D 2G)

Item No. 950100366

drive

DeVilbiss
HEALTHCARE



Flexibly using the bus and train for longer journeys? Not a problem with the Tourer. Thanks to the public transport authorization and the placard included with the scooter, it can be taken along on public transport on the spur of the moment. The ideal attendant and companion for city and countryside.

- The latest lighting technology reduces battery consumption
- Approved for transportation in designated buses in the public transport system with pictogram attached
- With the high-low mode, the speed can quickly be reduced to half the maximum speed by pressing a single button
- Identical model also available as 12 km/h scooter
- Lockable steering column which can be adjusted to suit
- Sleek, agile design

Technical Data

Order number	950100366
Colour	red
Overall width	64 cm
Overall length	120 cm
Overall height	115.5 cm
Seat width	46.5 cm
Seat depth	44 - 49 cm
Seat height	- cm
Width between armrest	43 - 55 cm
Speed	6 km/h

Battery capacity 2x	50 Ah
Motor capacity	500 W
Continuous-power	500 W
Gradeability	17 %
Clearance	8 cm
Front tires	11"
Rear tires	11"
Weight without battery	70,5 kg
Max. Load	136 kg

Accessories



Wheelchair / scooter ramp

Item No.STDS1097



Cane holder

Item No.950000028



Rain poncho

Item No.RT-SC001



Rain poncho

Item No.RT-SC004



Scooter garage, folding

Item No.SCGARAGE



Scooter bag for the backrest

Item No.RT-5026



Scooter saddle bag

Item No.RT-5022



Scooter accessory package

Item No.RT-SCPACK



Walking aid/rollator holder

Item No.950000089



Cane holder

Item No.950000029



Scooter garage

Item No.RT-4444