

SPARK (ST5D)

Item No. 950100420

drive

DeVilbiss
HEALTHCARE



With the Spark, your life becomes an experience. A number of intelligent functions and safety features ensure a perfect journey. Define your next journey yourself. The scooter is not only extremely safe but – with its snazzy design and the 8-inch aluminum rims with 13-inch low-profile tires – it's also a real eye-catcher.

- Powerful motor
- Individual mounting on the front wheels with adjustable degree of suspension for a ride to suit you
- With just one switch you can activate extra-slow mode and the maximum speed is reduced by half
- The volume of the beeps can be easily adjusted using the display
- Elegant 8" aluminium rims with 13" low profile tyres complete the timeless design
- The standard basket offers plenty of storage space
- The tiller is infinitely adjustable and can be adapted to the driver as required

Technical Data

Order number	950100420
Colour	anthracite
Overall width	68 cm
Overall length	144 cm
Overall height	136 cm
Seat width	51 cm
Seat depth	47 cm
Seat height	45 / 47.5 - 50 cm
Width between armrest	46 - 62 cm
Speed	15 km/h

Battery capacity 2x	80 Ah
Motor capacity	850 W
Continuous-power	850 W
Peak power	1300 W
Gradeability	10° / 17 %
Clearance	10 cm
Front tires	335 mm
Rear tires	335 mm
Weight without battery	138 kg
Max. Load	160 kg

Accessories



Wheelchair / scooter ramp

Item No.STDS1097



Rain poncho

Item No.RT-SC001



Rain poncho

Item No.RT-SC004



Scooter garage, folding

Item No.SCGARAGE



Scooter bag for the backrest

Item No.RT-5026



Scooter saddle bag

Item No.RT-5022



Scooter accessory package

Item No.RT-SCPACK



Walking aid/rollator holder

Item No.950000089



Cane holder

Item No.950000029



Scooter garage

Item No.RT-4444