

BLAZE (ST5D Plus)

Item No. 950100470

drive

DeVilbiss
HEALTHCARE



**Please note the changed item codes of our scooter portfolio for the international market. These change as follows:
Item No. 950100475-MD**

If you have any further questions, please contact our customer service.

Dynamic, sporty & sophisticated. The new Blaze is a crossover which was developed to meet your requirements. The innovative design puts the focus on the person. It is impressive not only because of the sophisticated details, but also because of its distinctive look. The 13-inch tires on stylish alloy rims complement the elegant design.

- Attentive driving with daytime running light as standard
- New, integrated LED lighting technology reduces power consumption
- Powerful 850 W motor
- 13" tyres on sleek design light metal aluminium rims ensure a comfortable drive
- The standard basket offers plenty of space to put your shopping, for example
- Non-slip throttles give a safe driving experience
- Volume of the beeps can be easily adjusted using the control panel
- Full suspension absorbs small bumps
- Extra-slow mode enables rapid, simple speed reduction to half the maximum speed

Technical Data

BLAZE (ST5D Plus)

Item No. 950100470



Order number	950100470
Colour	blue
Overall width	68 cm
Overall length	140 cm
Overall height	132.5 cm
Seat width	51 cm
Seat depth	47 cm
Seat height	- cm
Width between armrest	62 cm
Speed	15 km/h

Battery capacity 2x	80 Ah
Motor capacity	850 W
Gradeability	10° / 17 %
Clearance	9 cm
Front tires	335 mm
Rear tires	335 mm
Total weight	124 kg
Weight without battery	124 kg
Max. Load	160 kg

Accessories



Wheelchair / scooter ramp

Item No.STDS1097



Rain poncho

Item No.RT-SC001



Rain poncho

Item No.RT-SC004



Scooter garage, folding

Item No.SCGARAGE



Scooter bag for the backrest

Item No.RT-5026



Scooter saddle bag

Item No.RT-5022



Scooter accessory package

Item No.RT-SCPACK



Walking aid/rollator holder

Item No.950000089



Cane holder

Item No.950000029



Scooter garage

Item No.RT-4444