

SPARK 3 (ST5D 3 Wheel)

Item No. 950100503

drive

DeVilbiss
HEALTHCARE



Flexibly using the bus and train for longer journeys? Not a problem with the Tourer. Thanks to the public transport authorization and the placard included with the scooter, it can be taken along on public transport on the spur of the moment. The ideal attendant and companion for city and countryside.

- The latest lighting technology reduces battery consumption
- Approved for transportation in designated buses in the public transport system with pictogram attached
- With the high-low mode, the speed can quickly be reduced to half the maximum speed by pressing a single button
- Identical model also available as 12 km/h scooter
- Lockable steering column which can be adjusted to suit
- Sleek, agile design

Technical Data

Order number	950100503
Colour	anthracite
Overall width	62.5 cm
Overall length	143 cm
Overall height	116 cm
Seat height	- cm
Speed	15 km/h
Battery capacity 2x	80 Ah

Motor capacity	800 W
Gradeability	17 %
Clearance	10 cm
Front tires	13"
Rear tires	13"
Weight	113.7 kg
Max. Load	160 kg

SPARK 3 (ST5D 3 Wheel)

Item No. 950100503



Accessories



Wheelchair / scooter ramp

Item No.STDS1097



Cane holder

Item No.950000028



Rain poncho

Item No.RT-SC001



Rain poncho

Item No.RT-SC004



Scooter garage, folding

Item No.SCGARAGE



Scooter bag for the backrest

Item No.RT-5026



Scooter saddle bag

Item No.RT-5022



Scooter accessory package

Item No.RT-SCPACK



Walking aid/rollator holder

Item No.950000089



Cane holder

Item No.950000029



Scooter garage

Item No.RT-4444