

Scooter ST6

Item No. 950100611



- Simple and intuitive operation of the multifunction display
- Charger (12 Ah) with readout function
- Powerful 800 W motor
- Volume adjustment of the acoustic signals via key combination on the display
- Individually adjustable deluxe captain's seat
- Full suspension of all wheels for optimal ride comfort and safety when driving over obstacles
- Supplied as standard with bumper

Technical Data

Order number	950100611
Colour	red
Overall width	72 cm
Overall length	158 cm
Overall height	136 cm
Seat width	58 cm
Seat depth	56 cm
Seat height	- cm
Width between armrest	55 cm
Speed	15 km/h
Battery capacity 2x	80 Ah

Motor capacity	800 W
Continuous-power	800 W
Peak power	1600 W
Gradeability	9 %
Front tires	3.5 - 8
Rear tires	3.5 - 8
Total weight	171 kg
Weight without battery	120 kg
Max. Load	180 kg
PU pc.	1 Stk.

Scooter ST6

Item No. 950100611



Accessories



Rain poncho

Item No.RT-SC001



Rain poncho

Item No.RT-SC004



Scooter garage, folding

Item No.SCGARAGE



Scooter bag for the backrest

Item No.RT-5026



Scooter saddle bag

Item No.RT-5022



Walking aid/rollator holder

Item No.950000089



Cane holder

Item No.950000029



Scooter garage

Item No.RT-4444