

BLUSTER (ST6) 18 km/h

Item No. 950100700

drive

DeVilbiss
HEALTHCARE



Only for the international market.

Maximum power and top safety open up new perspectives for you! From the very first sight, Bluster meets all requirements. Its safe and convenient features ensure that you feel completely comfortable and can focus entirely on the enjoyment of driving during every outing.

- Well protected if it rains or the ground is muddy thanks to the splash guard
- The bumper provides protection and safety for the scooter
- High-quality 15" drive wheels with full suspension ensure the very best drive comfort
- The beep volume can be adjusted using a key sequence on the display
- Full cover of the rear wheel with integrated luminous backlights
- Sporty design and aerodynamic construction
- Scooter slows speed automatically when going round bends
- Powerful motor

Technical Data

Order number	950100700
Colour	anthracite
Overall width	72 cm
Overall length	160 cm
Overall height	136 cm
Seat width	58 cm
Seat depth	56 cm
Seat height	48 / 50.5 - 53 cm
Width between armrest	55 cm
Speed	18 km/h
Battery capacity 2x	80 Ah

Motor capacity	1000 W
Continuous-power	1000 W
Peak power	1600 W
Gradeability	16 %
Clearance	11 cm
Front tires	3.5 - 8
Rear tires	3.5 - 8
Weight without battery	120 kg
Max. Load	180 kg
PU Pal.	1 Pal.
PU pc.	1 Stk.

BLUSTER (ST6) 18 km/h

Item No. 950100700



Accessories



Cane holder

Item No.950000028



Rain poncho

Item No.RT-SC001



Rain poncho

Item No.RT-SC004



Scooter garage, folding

Item No.SCGARAGE



Scooter bag for the backrest

Item No.RT-5026



Scooter saddle bag

Item No.RT-5022




Walking aid/rollator holder

Item No.950000089



Cane holder

Item No.950000029



Top case, large

Item No.950100101



Scooter garage

Item No.RT-4444