

# ENVOY (pneumatic tires)

Item No. 950200100-50

**drive**

**DeVilbiss**  
HEALTHCARE



There is a good reason the Envoy is among the most popular scooters – its unique agility is impressive and guarantees plenty of driving pleasure. With numerous equipment features and a modern design, it guarantees mobility.

- With two rear view mirrors as standard for a good vision in town or in the countryside
- Very robust and strong with a maximum load of 160 kg
- Travel in public transport is possible with the pictogram delivered with or fitted on the product (Germany only)
- First-class, comfortable full suspension, as well as adjustment of all wheels to suit your needs
- Large foot area for greater freedom of movement when driving
- Chic, modern design
- Also available with pneumatic tyres for a comfortable ride over small obstacles and bumps

## Technical Data

|                       |                   |
|-----------------------|-------------------|
| Order number          | 950200100-50      |
| Colour                | gray metallic     |
| Overall width         | 58 cm             |
| Overall length        | 119 cm            |
| Overall height        | 117 - 122 cm      |
| Seat width            | 48 cm             |
| Seat depth            | 45 cm             |
| Seat height           | 46 / 48.5 - 51 cm |
| Width between armrest | 47 - 67 cm        |
| Speed                 | 6 km/h            |
| Battery capacity 2x   | 34/50 Ah          |


|                        |              |
|------------------------|--------------|
| Motor capacity         | 400 W        |
| Continuous-power       | 400 W        |
| Peak power             | 500 W        |
| Gradeability           | 9° / 16 %    |
| Clearance              | 8 cm         |
| Front tires            | 10" (260 mm) |
| Rear tires             | 10" (260 mm) |
| Weight without battery | 76 kg        |
| Max. Load              | 160 kg       |
| PU Pal.                | 1 Pal.       |

## Accessories



**Lap belt XL**

Item No.900300321



**Wheelchair / scooter ramp**

Item No.STDS1097



**Cane holder**

Item No.950000028



**Rain poncho**

Item No.RT-SC001



**Rain poncho**

Item No.RT-SC004



**Scooter garage L**

Item No.RT-4443



**Scooter garage, folding**

Item No.SCGARAGE



**Scooter bag for the backrest**

Item No.RT-5026



**Scooter saddle bag**

Item No.RT-5022



**Scooter accessory package**

Item No.RT-SCPACK



**Walking aid/rollator holder**

Item No.950000089



**Cane holder**

Item No.950000029



**Scooter garage**

Item No.RT-4444