

ENVOY (PU tires)

Item No. 950200150-36

drive

DeVilbiss
HEALTHCARE



- User weight up to 160 kg
- Full suspension for optimal ride comfort
- Modern design
- Simple and clear operation
- Rear-view mirrors on the left and right
- Supplied as standard with a front basket and cane holder
- Approved for carriage in intended local public transport buses with displayed pictogram
- With PU tires

Technical Data

Order number	950200150-36
Colour	gray metallic
Overall width	58 cm
Overall length	119 cm
Overall height	117 - 122 cm
Seat width	48 cm
Seat depth	45 cm
Seat height	46 / 48.5 - 51 cm
Width between armrest	47 - 67 cm
Speed	6 km/h
Battery capacity 2x	34/50 Ah


Motor capacity	400 W
Continuous-power	400 W
Peak power	500 W
Gradeability	9° / 16 %
Clearance	8 cm
Front tires	10" (260 mm)
Rear tires	10" (260 mm)
Weight without battery	76 kg
Max. Load	160 kg
PU Pal.	1 Pal.

Accessories



Lap belt XL

Item No.900300321



Wheelchair / scooter ramp

Item No.STDS1097



Cane holder

Item No.950000028



Rain poncho

Item No.RT-SC001



Rain poncho

Item No.RT-SC004



Scooter garage L

Item No.RT-4443



Scooter garage, folding

Item No.SCGARAGE



Scooter bag for the backrest

Item No.RT-5026



Scooter saddle bag

Item No.RT-5022



Scooter accessory package

Item No.RT-SCPACK



Walking aid/rollator holder

Item No.950000089



Cane holder

Item No.950000029



Scooter garage

Item No.RT-4444