

# Scooter BL350 Envoy

Item No. 950200150-50

**drive**

**DeVilbiss**  
HEALTHCARE



- User weight up to 160 kg
- Patented suspension for optimal driving comfort
- Modern design
- Easy and clear operation
- Mirrors left and right
- Extra-large wheels for maximum comfort and safe navigation of hazards

## Technical Data

|                       |               |
|-----------------------|---------------|
| Order number          | 950200150-50  |
| Colour                | gray metallic |
| Overall width         | 58 cm         |
| Overall length        | 119 cm        |
| Overall height        | 117 - 122 cm  |
| Seat width            | 48 cm         |
| Seat depth            | 45 cm         |
| Seat height           | - cm          |
| Width between armrest | 47 - 67 cm    |
| Speed                 | 6 km/h        |
| Battery capacity 2x   | 34/50 Ah      |


|                        |              |
|------------------------|--------------|
| Motor capacity         | 400 W        |
| Continuous-power       | 400 W        |
| Peak power             | 500 W        |
| Gradeability           | 9° / 16 %    |
| Clearance              | 8 cm         |
| Front tires            | 10" (260 mm) |
| Rear tires             | 10" (260 mm) |
| Weight without battery | 76 kg        |
| Max. Load              | 160 kg       |
| PU Pal.                | 1 Pal.       |

## Accessories



**Lap belt XL**

Item No.900300321



**Wheelchair / scooter ramp**

Item No.STDS1097



**Cane holder**

Item No.950000028



**Rain poncho**

Item No.RT-SC001



**Rain poncho**

Item No.RT-SC004



**Scooter garage L**

Item No.RT-4443



**Scooter garage, folding**

Item No.SCGARAGE



**Scooter bag for the backrest**

Item No.RT-5026



**Scooter saddle bag**

Item No.RT-5022



**Scooter accessory package**

Item No.RT-SCPACK



**Walking aid/rollator holder**

Item No.950000089



**Cane holder**

Item No.950000029



**Scooter garage**

Item No.RT-4444