

ENVOY (pneumatic tires)

Item No. 950200200-50

drive

DeVilbiss
HEALTHCARE



There is a good reason the Envoy is among the most popular scooters – its unique agility is impressive and guarantees plenty of driving pleasure. With numerous equipment features and a modern design, it guarantees mobility.

- With two rear view mirrors as standard for a good vision in town or in the countryside
- Very robust and strong with a maximum load of 160 kg
- Travel in public transport is possible with the pictogram delivered with or fitted on the product (Germany only)
- First-class, comfortable full suspension, as well as adjustment of all wheels to suit your needs
- Large foot area for greater freedom of movement when driving
- Chic, modern design
- Also available with pneumatic tyres for a comfortable ride over small obstacles and bumps

Technical Data

Order number	950200200-50
Colour	cherry-red
Overall width	58 cm
Overall length	119 cm
Overall height	117 - 122 cm
Seat width	48 cm
Seat depth	45 cm
Seat height	46 / 48.5 - 51 cm
Width between armrest	47 - 67 cm
Speed	6 km/h
Battery capacity 2x	34/50 Ah


Motor capacity	400 W
Continuous-power	400 W
Peak power	500 W
Gradeability	9° / 16 %
Clearance	8 cm
Front tires	10" (260 mm)
Rear tires	10" (260 mm)
Weight without battery	76 kg
Max. Load	160 kg
PU Pal.	1 Pal.

Accessories



Lap belt XL

Item No.900300321



Wheelchair / scooter ramp

Item No.STDS1097



Cane holder

Item No.950000028



Rain poncho

Item No.RT-SC001



Rain poncho

Item No.RT-SC004



Scooter garage L

Item No.RT-4443



Scooter garage, folding

Item No.SCGARAGE



Scooter bag for the backrest

Item No.RT-5026



Scooter saddle bag

Item No.RT-5022



Scooter accessory package

Item No.RT-SCPACK



Walking aid/rollator holder

Item No.950000089



Cane holder

Item No.950000029



Scooter garage

Item No.RT-4444