

FLAME 3 (HS925)

Item No. 950200310

drive

DeVilbiss
HEALTHCARE



- Comfortable seat
- Spacious storage compartment
- Adjustable steering column
- 15" large wheels
- Front and rear suspension
- LCD display
- LED lighting
- Comes standard with 2 mirrors

Technical Data

Order number	950200310
Overall width	69 cm
Overall length	159 cm
Overall height	127.5 cm
Seat width	45.6 cm
Seat height	- cm
Speed	15 km/h
Battery capacity 2x	50 Ah

Motor capacity	700 W
Gradeability	10 %
Front tires	370
Rear tires	370
Battery Voltage	12 V
Weight battery	30,5 kg
Weight without battery	69,5 kg
Max. Load	136 kg