



- Intelligent control panel with joystick
- Elegant black frame
- Ultralight and compact - perfect for travel and trips
- USB port for charging
- Easy-to-reach under-seat storage
- Removable and easy to clean seat
- Ventilated backrest
- Flip-up armrests
- Puncture-proof tires (rear wheels 8.5 inches)
- Suitable for indoor and outdoor use
- Can be loaded up to 115 kg
- With black transport bag including trolley function

Technical Data

Order number	ACFPC17BLK-EU
Colour	anthracite
Overall width	57 cm
Overall length	93 cm
Overall height	93 cm
Seat width	44 cm
Seat depth	42 cm
Seat height	- 49 cm
Speed	6 km/h
Battery Capacity	20 Ah
Continuous-power	250 W
Gradeability	10 %
Clearance	6,5 cm

Front tires	17.8 x 4.5 cm
Rear tires	21.6 x 5.2 cm
Folded width	57 cm
Folded height	32 cm
Folded length	79 cm
Battery Voltage	24 V
Total weight	18.5 kg
Weight battery	2,7 kg
Weight without battery	15,8 kg
Max. Load	115 kg
PU Pal.	6 Pal.
PU pc.	1 Stk.

Accessories



Travel battery Instafold/Airfold

Item No.950700004