



- Intelligent control panel with joystick
- Elegant black frame
- Ultralight and compact - perfect for travel and trips
- USB port for charging
- Easy-to-reach under-seat storage
- Removable and easy to clean seat
- Ventilated backrest
- Adjustable and flip-up armrests
- Puncture-proof tires (rear wheels 8.5 inches)
- Suitable for indoor and outdoor use
- Can be loaded up to 115 kg
- With black transport bag including trolley function

Technical Data

| | |
|------------------|---------------|
| Order number | ACFPC17BLK-EU |
| Colour | anthracite |
| Overall width | 57 cm |
| Overall length | 93 cm |
| Overall height | 93 cm |
| Seat width | 44 cm |
| Seat depth | 42 cm |
| Seat height | - 49 cm |
| Speed | 6 km/h |
| Battery Capacity | 20 Ah |
| Continuous-power | 250 W |
| Gradeability | 10 % |
| Clearance | 6,5 cm |

| | |
|------------------------|---------------|
| Front tires | 17.8 x 4.5 cm |
| Rear tires | 21.6 x 5.2 cm |
| Folded width | 57 cm |
| Folded height | 32 cm |
| Folded length | 79 cm |
| Battery Voltage | 24 V |
| Total weight | 18.5 kg |
| Weight battery | 2,7 kg |
| Weight without battery | 15,8 kg |
| Max. Load | 115 kg |
| PU Pal. | 6 Pal. |
| PU pc. | 1 Stk. |

Accessories



Travel battery Instafold/Airfold

Item No.950700004