

Hurricane Freedom edition, Purple

Item No. HCANE-PR

drive

DeVilbiss
HEALTHCARE



Our Hurricane walking stick folds and unfolds in seconds. Additionally it is not only a secure support, but also particularly chic and available in various fresh colours.

- Available in 4 different colours
- 8 convenient height adjustments
- Comfortable vinyl handgrip with wrist strap
- Non-slip and large-surface rubber foot
- Self-standing cane thanks to his large rubber foot

Technical Data

Order number	HCANE-PR
Handle height	77 - 95 cm
Total weight	500 g

Max. Load	159 kg
PU pc.	25 Stk.