

Fingertip pulse oximeter HbO-Smart

Item No. HbO-Smart



The new HbO-Smart device is the latest fingertip pulse oximeter from Drive DeVilbiss intended for use in homes or hospitals for non-invasive measurement of oxygen saturation, pulse rate, and perfusion index. With settings for user alerts and four optional screen displays, this new, elegant device offers everything you would expect from a modern oximeter that can be used for both children and adults

- OLED display with heart rate bar and plethysmography waveform
- Four different screen settings
- Equipped with a low battery alarm and automatic shutoff after 8 seconds if no signal is detected
- Low power consumption (2 AAA batteries for 25 hours of operation)

Technical Data

Order number	HbO-Smart
SP02 measurement	70 - 99%
Overall width	3.5 cm
Overall length	6.5 cm
Overall height	3.0 cm

Pulse rate measurement	30 - 240 BPM
Pulsfrequenz Genauigkeit	80% - 99% ± 2%
Total weight	6 g
Weight	6 g