




- Intelligent control panel with joystick
- Chassis in metallic silver
- Rims in silver-black
- Removable padded backrest and seat cushion
- Front wheel suspension
- Standard pocket under the seat
- Height adjustable footfalls
- Foldable armrests
- lap belt
- Folded in a few seconds

Technical Data

Order number	IFPC17SIL-FR
HMV-Nr.	18.50.07.0011
Colour	white
Overall width	62 cm
Overall length	105 cm
Overall height	99 cm
Seat width	44 cm
Seat depth	45 cm
Seat height	- 55 cm
Speed	6 km/h
Battery Capacity	20 Ah
Motor capacity	200 W
Gradeability	10 %

Clearance	7,5 cm
Front tires	20 x 5 cm
Rear tires	30.5 x 5.7 cm
Folded width	62 cm
Folded height	40 cm
Folded length	79 cm
Total weight	27 kg
Weight battery	2,7 kg
Weight without battery	24,3 kg
Max. Load	120 kg
PU Pal.	4 Pal.
PU pc.	1 Stk.

Accessories



Travel battery Instafold/Airfold

Item No.950700004