

# Hikoneb 908 DC

Item No. M000034



The reliable Hikoneb 908 DC is a high performance, large volume ultrasonic nebuliser suitable for room humidification, sputum induction and highly effective aerosol therapy.

Designed with the patient and end user in mind, the Hikoneb 908 device is simple to use, and ideal for use in both clinical and care home environments.

- Built in medication timer up to 90 Min.
- Air-flow output with 4 adjustable settings
- Adjustable nebulisation rate
- Digital control panel, allows precise and easy operation
- Heated tube option and 5-feet trolley included
- Low water level audible alarm

Fully equipped with a heated tube, a 5-feet trolley, medical cups and all accessories, the Hikoneb 908 DC is ready for use, right out of the box.

## Technical Data

Order number	M000034
Overall width	22.6 cm
Overall length	10.8 cm
Overall height	29 cm
Air flow max.	19 l/min
Nebulizer rate	2,3 ml/min

Particle size MMAD	< 4 µm
Power Supply	50 Hz
Ultrasonic frequency	1.7 Mhz
Weight	3 kg
PU pc.	1 Stk.