

Compressor nebulizer AirForce One

Item No. MQ5890-EU



- Compact compressor nebulizer for reliable and inexpensive aerosol therapy
- Simple to set up and operate by the user
- Reliable and inexpensive therapy device
- Suitable for daily use
- Quiet operation
- 2-year warranty
- Integrated carrying handle for additional comfort
- Extensive accessories and storage bag included in delivery
- Equally suitable for children and adults

Technical Data

Order number	MQ5890-EU
Colour	white
Overall width	14 cm
Overall length	19 cm
Overall height	9 cm
Aerosol output	0,45 ml
Air flow max.	5.1 l/min
Nebulizer rate	0,3 ml/min

Particle size MMAD	< 5 µm
Power Supply	50 Hz
Sound level	52 dBa
Weight	1.5 kg
Weight without accessories	1.5 kg
PU Pal.	168 Pal.
PU pc.	4 Stk.