



- Integrated battery with LED status indicators
- Light as a feather
- Particularly pleasant tactile qualities
- Simple handling and cleaning
- Low operating noise
- One of the smallest inhalation devices on the market

The AirForce Mini is small, light, and handy, making it an ideal companion for effective inhalation therapy at home and on the go.

The compressor technology that is used produces a very high percentage of respirable aerosols, making the AirForce Mini outstandingly suitable for treating diseases of the lower and upper airways. The AirForce Mini has been designed for use both by children and adolescents and by adults.

Included in delivery:

- AirForce Mini compressor
- Mouthpiece with exhalation valve
- Air tube
- Replacement filters (5 pieces)
- Adult mask
- Pediatric mask
- Nosepiece
- Micro USB cable
- Nozzle nebulizer
- Instructions for use
- Carrying bag and power supply unit with plug

Technical Data

Order number	MQ6010
Overall width	6.8 cm
Overall length	11,1 cm
Overall height	4.5 cm
Aerosol output	0,15 ml
Air flow max.	0.15 l/min

Nebulizer rate	0,25 ml/min
Particle size MMAD	4,12 µm
Power Supply	50 Hz
Weight	217 g
Weight without accessories	217 g
PU pc.	6 Stk.