



With the Raven, you are free and mobile at all times when selecting your activities. This is because the maneuverable scooter's modern design is not just attractive but also allows it to be driven very easily.

- Can be disassembled quickly without the need for tools and is easy to store
- Low net weight
- Clear control panel guarantees easy operation
- Small turning circle ensures exceptional manoeuvrability
- Standard pneumatic tyres absorb bumps
- Anti-tip wheels ensure safe driving over small curbs or landings
- Seat height can be adjusted to 3 heights (54 cm, 56.5 cm, 59 cm)

Technical Data

Order number	MS025RD-G-AIR
Colour	red
Overall width	48 cm
Overall length	108 cm
Overall height	91 cm
Seat width	43 cm
Seat depth	36 cm
Seat height	44 / 56.5 - 59 cm
Width between armrest	40 - 50 cm
Speed	6 km/h
Battery capacity 2x	20 Ah

Motor capacity	270 W
Continuous-power	270 W
Peak power	500 W
Gradeability	8° / 14 %
Clearance	6 cm
Front tires	200 x 50 mm
Rear tires	200 x 50 mm
Weight without battery	34 kg
Max. Load	130 kg
PU Pal.	1 Pal.

Accessories



Rain poncho

Item No. RT-SC001



Rain poncho

Item No. RT-SC004



Scooter garage L

Item No. RT-4443



Scooter garage, folding

Item No. SCGARAGE



Scooter bag for the backrest

Item No. RT-5026



Scooter saddle bag

Item No. RT-5022



Scooter accessory package

Item No. RT-SCPACK



Scooter garage M

Item No. RT-2993



Cane holder

Item No. 950000029