

Pneumatic oxygen conserver Bonsai

Item No. OM-824

drive

DeVilbiss
HEALTHCARE



- Potential saving ratio of up to 6:1
- 7 settings in pulsed mode
- 2 l/min with continuous flow
- Maintenance-free
- Light and compact
- Battery-free operation
- Inexpensive and efficient
- High respiratory frequency range of up to 40 BPM with the same bolus volume to reduce the risk of desaturation during patient activity

Technical Data

Order number	OM-824
Colour	black
Overall width	6.5 cm

Overall height	6.5 cm
Connector type	G3/4. DIN ISO 477 Nr. 9
Total weight	300 g

Pneumatic oxygen conserver Bonsai

Item No. OM-824



Accessories



Trolley for oxygen bottle

Item No.210000001



Bag for oxygen bottle

Item No.EX3000D-651