

Electronic oxygen conserver Evolution

Item No. OM-924



- Potential saving ratio of up to 5:1
- 7 settings in demand mode
- 2 l/min with continuous flow
- Light and compact
- Inexpensive and efficient
- High respiratory frequency range of up to 40 BPM with the same bolus volume to reduce the risk of desaturation during patient activity
- Responds to both weak and strong inhalation, making it suitable for a broad range of patients
- Automatic shutoff when not in use, thus extending the battery life
- Evolution Motion with integrated motion sensor for detecting patient activity
- Long battery life

Technical Data

Order number	OM-924
Colour	black
Overall width	7.9 cm

Overall height	6.4 cm
Connector type	G3/4. DIN ISO 477 Nr. 9
Total weight	600 g

Electronic oxygen conserver Evolution

Item No. OM-924

drive

DeVilbiss
HEALTHCARE

Accessories



Trolley for oxygen bottle

Item No.210000001



Bag for oxygen bottle

Item No.EX3000D-651