

Disposable oxygen conserver Oxymizer

Item No. P-224



- Can be used with a continuous flow of up to 15 l/min.
- Comfortable alternative to an oxygen mask
- Facilitates the delivery of continuous high-flow oxygen therapy in a home care, hospital, hospice or long-term care setting
- Costs lower thanks to reduced oxygen use
- Lower flow rates reduce the necessity of a humidifier
- Enables the patient to eat, drink, and talk
- Compatible with compressed gas cylinders, oxygen concentrators, and liquid oxygen systems
- Disposable, does not require batteries
- Clinically proven efficiency when resting, sleeping, and during physical activity
- Oxygen saving ratio up to 4:1

Technical Data

Order number

P-224

Weight

8.5 g