

# Wheelchair pushing aid PowerStroll

Item No. PWCPP011



- Puncture-proof drive wheels
- Can be installed in just a few steps
- Suitable for almost all standard wheelchairs
- The PowerStroll is installed in a secure position beneath the wheelchair
- The charging time is eight hours with no risk of overcharging
- Clear and easy-to-use control
- Continuously variable speed regulation up to 6 kph forwards and backwards
- Suitable for wheelchairs with a seat width of 37 to 58 cm
- Including telescopic bar (Item No. 960000003) suitable for wheelchairs with a seat width of 37 to 58 cm

## Technical Data

Order number	PWCPP011
Colour	black
Overall width	25 cm
Overall length	41 cm
Seat height	- cm
Speed	6 km/h
Battery capacity 2x	12 Ah

Motor capacity	150 W
Gradeability	10.5 %
Total weight	19.9 kg
Weight battery	9,7 kg
Weight without battery	10,2 kg
Max. Load	135 kg