

Oxygen meter UltraMaxO₂ analyzer

Item No. R221P11-001

drive

DeVilbiss
HEALTHCARE



The UltraMaxO₂ analyzer was developed to check the oxygen concentration at a glance, including the flow and output pressure of oxygen concentrators. Its advanced design offers high performance and reliability. The UltraMaxO₂ oxygen analyzer offers the following advantageous characteristics:

- Precise measurement of the O₂ flow and pressure
- Integrated oxygen sensor; no replacement of the oxygen sensor necessary
- No calibration on site
- Self-diagnosis
- Shielded and reinforced sample gas inlet
- Including silicone sleeve for increased durability
- 36-month warranty
- Maintenance-free
- Manufactured in the USA

Technical Data

Order number	R221P11-001
Overall width	8.6 cm
Overall height	12.9 cm
Overall depth	2.5 cm

Flow rate	0 - 10 l/min
Total weight	181 g
Weight	181 g