

Patient lift Hermes

Item No. SR2310000



- The powerhouse
- Heavy-duty patient lifter, lifting capacity up to 250 kg
- Extra-reinforced steel frame
- Electric chassis spreading for effortless maneuvering
- Low base frame height for easy movement beneath furnishings
- Four-point transport frame with protective covering
- Manual emergency lowering
- Double, full-rubber casters, diameter 100 mm

[Click here to download our checklist for reprocessing and maintenance \(dealer login required\)](#)

Technical Data

Order number	SR2310000
Colour	gray
Overall width	74 cm
Overall length	136 cm
Expansion area	60 - 110 cm

Height difference	51 - 200 cm
Turning radius	300 cm
Clearance	11,4 cm
Total weight	57 kg
Max. Load	250 kg

Patient lift Hermes

Item No. SR2310000



Accessories



Patient scale for Hermes

Item No. SR2321000