

Patient lift Samsoft 175

Item No. SR2870000



- Versatile and powerful
- Can be broken down into two sections for easy transport and handling without tools
- Reinforced steel frame for maximum safety
- Large adjustment range of the chassis spread for optimal stability and easy avoidance of obstacles
- Chassis spread with lock
- Lightweight, compact, and folding
- Easy and quick battery charging
- Supplied with a sling as standard, optionally in size M or L

[Click here to download our checklist for reprocessing and maintenance \(dealer login required\)](#)

[Important, safety-relevant information about this product \(dealer login required\)](#)

Technical Data

Order number	SR2870000
Colour	gray
Overall width	63 cm
Overall length	173 cm
Overall height	136 cm
Chassis length closed	54
Expansion area	52 - 72 cm

Height difference	55 - 189 cm
Turning radius	144 cm
Clearance	9 cm
Total weight	38 kg
Max. Load	175 kg
PU Pal.	4 Pal.

Patient lift Samsoft 175

Item No. SR2870000

drive

DeVilbiss
HEALTHCARE

Accessories



Patient scale for Samsoft 175

Item No. SR2820400