

CROW (ST3)

Item No. ST3406100

drive

DeVilbiss
HEALTHCARE



Does planning your day freely and flexibly sound like a dream to you? The Crow makes it easy, because the maneuverable scooter not only features a modern design but also simple and flexible handling. The Crow is ideally suited for traveling and excursions because the scooter can be disassembled in just a few easy steps without tools and stowed quickly in the car.

- Can be disassembled without tools
- Small and maneuverable
- Lots of legroom despite the compact size
- Front LEDs for good visibility
- Rear suspension for a comfortable ride
- Available in the colors red and anthracite

Technical Data

Order number	ST3406100
Colour	red
Overall width	59 cm
Overall length	112 cm
Overall height	93.5 cm
Seat width	46 cm
Seat height	- cm
Speed	8 km/h
Battery capacity 2x	22 Ah

Motor capacity	250 W
Gradeability	14 %
Clearance	52 cm
Front tires	225 mm
Rear tires	225 mm
Total weight	55 kg
Weight	55 kg
Weight battery	15,5 kg
Max. Load	136 kg