

# CROW (ST3)

Item No. ST3406300

**drive**

**DeVilbiss**  
HEALTHCARE



Does planning your day freely and flexibly sound like a dream to you? The Crow makes it easy, because the maneuverable scooter not only features a modern design but also simple and flexible handling. The Crow is ideally suited for traveling and excursions because the scooter can be disassembled in just a few easy steps without tools and stowed quickly in the car.

- Can be disassembled without tools
- Small and maneuverable
- Lots of legroom despite the compact size
- Front LEDs for good visibility
- Rear suspension for a comfortable ride
- Available in the colors red and anthracite

## Technical Data

Order number	ST3406300
Colour	anthracite
Overall width	59 cm
Overall length	112 cm
Overall height	93.5 cm
Seat width	46 cm
Seat height	- cm
Speed	8 km/h
Battery capacity 2x	22 Ah
Motor capacity	250 W

Continuous-power	250 W
Gradeability	14 %
Clearance	52 cm
Front tires	225 mm
Rear tires	225 mm
Total weight	55 kg
Weight	55 kg
Weight battery	15,5 kg
Max. Load	136 kg