

Transport chair Travelite

Item No. TC005

drive

DeVilbiss
HEALTHCARE



- The ideal travel companion for car, plane, and train
- Stable, comfortable, and fits through any door
- Immediately ready to use and accompanies you wherever you want
- The soft seat surface and backrest offer comfort
- Including carrying bag for easy transport
- With a total weight of just 9.3 kg, ideal as hand baggage
- Folded in just a few easy steps
- Ensures maximum mobility

Technical Data

Order number	TC005
Colour	black-silver
Overall width	57 cm
Overall length	99 cm
Overall height	94 cm
Seat width	46 cm
Seat depth	43 cm

Seat height	- 48 cm
Folded height	29 cm
Folded width	27 cm
Folded length	77 cm
Total weight	9.3 kg
Max. Load	115 kg
PU Pal.	24 Pal.

Transport chair Travelite

Item No. TC005



Accessories



Carrying bag with transport rollers

Item No.900600800