

Rollator Tri Walker

Item No. TW011R

drive

DeVilbiss
HEALTHCARE



- Can be used indoors and outdoors
- Easy to transport
- Ergonomic, height-adjustable handles
- Passive lighting at the rear and front
- Supplied as standard with basket and tray

Technical Data

Order number	TW011R
Colour	red
Overall width	61 cm
Overall length	60 cm
Handle height	84.5 - 94.5 cm
Seat height	- cm
Width between handles	37 cm
Width between rear wheels	51 cm

Folded width	25 cm
Folded height	82,5 cm
Folded length	65 cm
Total weight	6.9 kg
Weight	6.2 kg
Weight without accessories	6.2 kg
Max. Load	130 kg
PU Pal.	12 Pal.

Rollator Tri Walker

Item No. TW011R



Accessories



One-hand brake

Item No.710000081



Basket

Item No.710500007



Tray

Item No.710500006



Bag

Item No.710500005