

Quad cane L

Item No. WS016L



- The quad cane offers maximum safety and ensures optimal stability
- The ergonomic hand grip sits comfortably and securely in the hand
- The height can be individually adjusted to 11 positions to suit the user
- The grip thrills with its high-quality finish

from week 21/22: Please note the change of article number to WS026L

Technical Data

Order number	WS016L
Colour	silver
Handle height	70 - 95 cm
Total weight	1.3 kg

Weight	1.3 kg
Max. Load	115 kg
PU pc.	4 Stk.