



Patented, detachable backrest

Easy to remove backrest

Easy care seat with practical hygiene cutout

Ergonomic backrest with inclination up to 50° degrees

Twice as safe: more stability thanks to two racks

Stable: patented 3D suction cups

Find more info on our website



Total width	70 cm
Seat width	38 cm
Seat depth	52,5 cm
Seat height	6,5 - 48 cm
Seat plate total width	64 cm
Backrest height	60 cm
Backrest inclination angle	10° - 50°
Base plate (L x B)	60 x 30 cm
Weight	9,5 kg
Max. load	140 kg

## YOU CAN RELY ON DRIVE DEVILBISS

We are the pioneers of the battery-powered bath lift and for over 30 years we have been a reliable partner for our customers.

See for yourself!



YOUR DEALER:

Drive Medical GmbH & Co. KG • Leutkircher Straße 44 • 88316 Isny/Allgäu  
 T +49 (0) 7562 9724 0 • info@drivemedical.de • www.drivedevilbiss.de

Stand: 12.08.2022  
 Irrtum und Änderungen vorbehalten



# BATHING WITHOUT LIMITS



BELLAVITA  
 NOVA

From Allgäu

... INTO THE WORLD

### TRADITION

As inventors of the first battery-powered bath lift, we can look back on more than 30 years of experience. And that you can see and feel in every detail of our **BELLAVITA NOVA**.

We are particularly proud of the production and development of all our Bellavita models at our location in the Allgäu in South-Germany. We attach great importance to working with regional suppliers.

Our bath lifter is a sustainable quality product „Made in Germany“ that we deliver from Germany to the whole world.

### INNOVATION

The **BELLAVITA NOVA** combines excellent functionality with technical finesse and an ergonomically handy design.

All of the experience from 30 years of bath lift development has been combined to create an innovative product that is extremely practical.

... INTO YOUR BATHROOM

#### Lightweight and portable:

Thanks to its low weight, the detachable **BELLAVITA NOVA** can be quickly placed in standard bathtubs.

#### Comfortable and practical:

Thanks to its maximum seat height of 48 cm, no additional adapter is required. The side flaps allow for easy transfer. Optionally, these are also available in wider versions.

#### Twice as safe:

The two patented toothed racks guarantee safe raising and lowering. Unique 3D suction cups ensure a robust and stable stand.

#### Infinitely variable fun:

The **BELLAVITA NOVA** can be smoothly lowered and raised using the buoyant hand control.

# BATHING WITHOUT LIMITS

... is possible again thanks to the **BELLAVITA NOVA**. People who have to do without a soothing bath due to physical limitations or various types of complaints can experience relaxation for their muscles and mind again. This is exactly why we developed the bath lift **BELLAVITA NOVA**.

EXPERIENCE THAT PAYS OFF FOR YOU

With the **BELLAVITA NOVA** experience meets innovation. The result: a bath lift that has not existed before. Equipped with a variety of practical details and patented innovations, it shines with benefits that guarantee unlimited bathing.

Many heads and hands were involved in the development of the **BELLAVITA NOVA**. But the heart of the project beats in our technician Robert Kempter. He was already instrumental in developing the very first battery-powered bath lift. His experience, our know-how and a lot of heart and soul you will find in our newest star in the bathing sky.



” ROBERT KEMPTER

„To get even better products - that is always my goal. I am sure that we succeeded with the new **BELLAVITA NOVA**“

“



Developed by

*R. Kempter*

Drive DeVilbiss Technician

